

POLISHING INSTRUMENTS
FOR THE DENTAL PRACTICE

ROTA
DIATECH





COMPREPOL ULTRA				
ISO Ø 1/10 mm	030	050	060	100
Head L in mm	7.5	10	7.5	2.5
Polisher type	2101RA	2103RA	2111RA	2127RA
RA 10 pcs	230055AA	230059AA	230063AA	230069AA
RA 50 pcs		230060AA	230064AA	

COMPOSHINE ULTRA				
ISO Ø 1/10 mm	030	050	060	100
Head L in mm	7.5	10	7.5	2.5
Polisher type	2201RA	2203RA	2211RA	2227RA
RA 10 pcs	230072AA	230076AA	230080AA	230086AA
RA 50 pcs		230077AA	230081AA	

Composite, Compomere ① ②



COMPREPOL PLUS				
ISO Ø 1/10 mm	030	040	060	100
Head L in mm	7.5	10	7.5	2.5
Polisher type	2301RA	2303RA	2311RA	2327RA
RA 5 pcs	230291AA	230293AA	230295AA	230297AA
RA 10 pcs		230294AA	230296AA	

COMPOSHINE PLUS				
ISO Ø 1/10 mm	030	040	060	100
Head L in mm	7.5	10	7.5	2.5
Polisher type	2401RA	2403RA	2411RA	2427RA
RA 5 pcs	230299AA	230301AA	230303AA	230305AA
RA 10 pcs		230302AA	230304AA	

Composite, Compomere ① ②



DIASHINE COMPOMANT PLUS				
ISO Ø 1/10 mm	030	040	060	100
Head L in mm	7.5	10	7.5	2.5
Polisher type	9101RA	9103RA	9111RA	9121RA
RA 5 pcs	230264AA	230265AA	230267AA	230271AA
RA 10 pcs	230415AA	230416AA		

UNISHINE			
ISO Ø 1/10 mm	055	070	090
Head L in mm	15	10	8.5
Polisher type	4203RA	4211RA	4216RA
RA 10 pcs	230203AA	230211AA	230215AA

Composite, Compomere ① ②



PREBROWN			
ISO Ø 1/10 mm	040	045	060
Head L in mm	7.5	10	9.5
Polisher type	1101RA	1103RA	1111RA
RA 10 pcs	230001AA	230005AA	230013AA
RA 50 pcs		230006AA	230014AA

HIGREEN			
ISO Ø 1/10 mm	040	045	060
Head L in mm	7.5	10	9.5
Polisher type	1201RA	1203RA	1211RA
RA 10 pcs	230028AA	230032AA	230040AA
RA 50 pcs			230041AA

Precious Metal Alloy, Amalgam ① ②



CERAPRESHINE		
ISO Ø 1/10 mm	050	070
Head L in mm	10	10
Polisher type	3201RA	3211RA
RA 10 pcs	230114AA	230122AA

CERASHINE		
ISO Ø 1/10 mm	050	070
Head L in mm	10	10
Polisher type	3301RA	3311RA
RA 10 pcs	230139AA	230147AA

ZIRCOPOL PLUS		
ISO Ø 1/10 mm	040	060
Head L in mm	10	7.5
Polisher type	9703RA	9711RA
RA 10 pcs	230379AA	230381AA

ZIRCOSHINE PLUS		
ISO Ø 1/10 mm	040	060
Head L in mm	10	7.5
Polisher type	9803RA	9811RA
RA 10 pcs	230391AA	230393AA

Ceramic ① → ②

CAD/CAM Restoration, Zirconia ① → ②




CERAFIN PLUS			
ISO Ø 1/10 mm	040	060	110
Head L in mm	10	7.5	2
Polisher type	9303RA	9311RA	9331RA
RA 5 pcs	230317AA	230319AA	230325AA
RA 10 pcs	230318AA		

CERAPRESHINE PLUS				
ISO Ø 1/10 mm	030	040	060	110
Head L in mm	7.5	10	7.5	2
Polisher type	9401RA	9403RA	9411RA	9431RA
RA 5 pcs	230331AA	230335AA	230337AA	230343AA

CERASHINE PLUS				
ISO Ø 1/10 mm	030	040	060	110
Head L in mm	7.5	10	7.5	2
Polisher type	9501RA	9503RA	9511RA	9531RA
RA 5 pcs	230349AA	230353AA	230355AA	230361AA
RA 10 pcs		230354AA		

Ceramic ① → ② → ③



BRUSH			
ISO Ø 1/10 mm	060	040	040
Head L in mm	5.5	6.5	6.5
Polisher type	7001RA	7002RA	7003RA
RA 5 pcs	232002AA	232004AA	232006AA
RA 10 pcs	232001AA	232003AA	232005AA

MANDREL	SWISSFLEX DISCS				
Length	Grit size	C/70µm	M/50µm	F/30µm	UF/5µm
25.4	Pack size	80 pcs	80 pcs	80 pcs	80 pcs
702RA	Polisher type	704UM ¹ /703UM ²	706UM ¹ /705UM ²	708UM ¹ /707UM ²	710UM ¹ /709UM ²
232018AA	Head Ø: 9 mm	232020AA ¹	232022AA ¹	232024AA ¹	232026AA ¹
	Head Ø: 13 mm	232019AA ²	232021AA ²	232023AA ²	232025AA ²

Brush Composite, Amalgam, Glass Ionomer Cement, Semi-precious Metal, Precious Metal



SWISSFLEX STRIPS	
Length x Width	175mm x 3.5mm
Pack size	80 pcs
700UM: C(70µm)/M(5µm)	232016AA
701UM: F(30µm)/UF(5µm)	232017AA

See SwissFlex Discs

DIAStrip	
Length x Width	80mm x 2.5mm
Polisher type	600UM
5 pcs	232014AA
10 pcs	232015AA




Order information

ISO number
(material, shank, shape, grit size,
variable diameter)

Head diameter in 1/10 mm
Head length
Polisher description
Packaging sizes available

ISO 658 204 243 523 ...



Illustration

ISO Ø 1/10 mm	050
Head L in mm	10
Polisher type	2103RA
RA 10 pcs	230059AA
RA 50 pcs	230060AA

Order numbers in respective
packaging size

Head Shapes

ISO Head Shape	Description	ISO Head Shape	Description
014	Small Wheel	243 040 10	Point
030	Cup	243 045-055 010	Flame
042	Wheel	273	Dome
044	Wheel	274	Bullet
107	Cylinder	292	Torpedo
114	Cylinder	303	Large Lens
173	Taper	304	Small Lens
243 030-040 7.5	Minipoint	372	Disc

Recommended Speed Ranges

Polishing System	RPM	Polishing System	RPM
COMPREPOL ULTRA	5'000-10'000	CERAPRESHINE	5'000-10'000
COMPOSHINE ULTRA	5'000-10'000	CERASHINE	5'000-10'000
COMPREPOL PLUS	3'000-8'000	ZIRCOPOL PLUS	5'000-10'000
COMPOSHINE PLUS	3'000-8'000	ZIRCOSHINE PLUS	5'000-10'000
DIASHINE COMPOMANT PLUS	3'000-8'000	CERAFIN PLUS	5'000-10'000
UNISHINE	5'000-10'000	CERAPRESHINE PLUS	5'000-10'000
PREBROWN	5'000-10'000	CERASHINE PLUS	5'000-10'000
HIGREEN	5'000-10'000	BRUSH	5'000-10'000
		SWISSFLEX	10'000-16'000

CH Coltene/Whaledent AG
Feldwiesenstr. 20
9450 Altstätten / Switzerland
T + 41 71 757 53 00
F + 41 71 757 53 01
info.ch@coltene.com

UK Coltene/Whaledent Ltd.
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF / UK
T +44 1444 235486
F +44 1444 870640
info.uk@coltene.com

US Coltene/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223 / USA
T + 1 330 916 8800
F + 1 330 916 7077
info.us@coltene.com

© 2016 Coltene/Whaledent AG – www.coltene.com



COLTENE

